



Filename: 2021-04-06 - 15.16.45.D5000



Default image (Contrast 100%)

Sample Info


Well	Conc. [ng/μl]	Sample Description	Alert	Observations
L1	39.6	Ladder	⚠	Caution! Expired ScreenTape device; Ladder
B1	1.66	9 (2550)	⚠	Caution! Expired ScreenTape device
D1	17.8	11 (2877)	⚠	Caution! Expired ScreenTape device
E1	12.9	12 (1732)	⚠	Caution! Expired ScreenTape device
L2	38.8	Ladder	⚠	Caution! Expired ScreenTape device; Ladder
F1	12.1	13 (1329)	⚠	Caution! Expired ScreenTape device
G1	23.1	14 (1563)	⚠	Caution! Expired ScreenTape device
A2	20.0	16 (1728)	⚠	Caution! Expired ScreenTape device
B2	28.5	17 (2012)	⚠	Caution! Expired ScreenTape device
C2	25.0	18 (1051)	⚠	Caution! Expired ScreenTape device

L1: Ladder**Sample Table**

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
L1	39.6	Ladder		Caution! Expired ScreenTape device; Ladder

Peak Table


Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	5.91	-	606	-		Lower Marker
100	3.76	-	57.9	9.50		
250	4.24	-	26.1	10.72		
400	4.35	-	16.7	11.00		
600	4.58	-	11.7	11.57		
1000	4.50	-	6.93	11.38		
1500	4.53	-	4.64	11.44		
2500	4.43	-	2.73	11.19		
3500	4.56	-	2.00	11.52		
5000	4.63	-	1.42	11.69		
10000	3.25	3.25	0.500	-		Upper Marker

B1: 9 (2550)**Sample Table**

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
B1	1.66	9 (2550)		Caution! Expired ScreenTape device

Peak Table


Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	6.45	-	661	-		Lower Marker
569	1.66	-	4.50	100.00		
10000	3.25	3.25	0.500	-		Upper Marker

D1: 11 (2877)**Sample Table**

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
D1	17.8	11 (2877)	⚠	Caution! Expired ScreenTape device

Peak Table


Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	7.01	-	719	-		Lower Marker
297	17.8	-	92.5	100.00		
10000	3.25	3.25	0.500	-		Upper Marker

E1: 12 (1732)**Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
E1	12.9	12 (1732)		Caution! Expired ScreenTape device

Peak Table


Size [bp]	Calibrated Conc. [ng/ul]	Assigned Conc. [ng/ul]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	6.19	-	635	-		Lower Marker
286	12.9	-	69.2	100.00		
10000	3.25	3.25	0.500	-		Upper Marker

L2: Ladder**Sample Table**

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
L2	38.8	Ladder		Caution! Expired ScreenTape device; Ladder

Peak Table


Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	6.04	-	619	-		Lower Marker
100	3.77	-	58.0	9.73		
250	4.17	-	25.7	10.76		
400	4.22	-	16.2	10.88		
600	4.45	-	11.4	11.48		
1000	4.41	-	6.78	11.37		
1500	4.42	-	4.53	11.41		
2500	4.37	-	2.69	11.27		
3500	4.55	-	2.00	11.74		
5000	4.41	-	1.36	11.37		
10000	3.25	3.25	0.500	-		Upper Marker

F1: 13 (1329)**Sample Table**

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
F1	12.1	13 (1329)	⚠	Caution! Expired ScreenTape device

Peak Table


Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	6.86	-	704	-		Lower Marker
359	12.1	-	51.8	100.00		
10000	3.25	3.25	0.500	-		Upper Marker

G1: 14 (1563)**Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
G1	23.1	14 (1563)		Caution! Expired ScreenTape device

Peak Table


Size [bp]	Calibrated Conc. [ng/ul]	Assigned Conc. [ng/ul]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	6.72	-	689	-		Lower Marker
315	23.1	-	113	100.00		
10000	3.25	3.25	0.500	-		Upper Marker

A2: 16 (1728)**Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
A2	20.0	16 (1728)		Caution! Expired ScreenTape device

Peak Table


Size [bp]	Calibrated Conc. [ng/ul]	Assigned Conc. [ng/ul]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	6.79	-	696	-		Lower Marker
387	20.0	-	79.7	100.00		
10000	3.25	3.25	0.500	-		Upper Marker

B2: 17 (2012)**Sample Table**

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
B2	28.5	17 (2012)		Caution! Expired ScreenTape device

Peak Table

Size [bp]	Calibrated Conc. [ng/ul]	Assigned Conc. [ng/ul]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	6.97	-	715	-		Lower Marker
353	28.5	-	124	100.00		
10000	3.25	3.25	0.500	-		Upper Marker

C2: 18 (1051)**Sample Table**

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
C2	25.0	18 (1051)	⚠	Caution! Expired ScreenTape device

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	8.10	-	831	-		Lower Marker
431	25.0	-	89.2	100.00		
10000	3.25	3.25	0.500	-		Upper Marker